


Series 100 Argon Gas Chromatograph

AGC
INSTRUMENTS

Gas Chromatography since 1965



..... Analysis of Impurities in Argon to ppb levels

Highly recommended for the quality control of Argon Gas. We can measure Impurities of H₂, O₂, N₂, CH₄, CO & CO₂ to ppb levels with ease. Easy to use with minimal training required, ensure the Series 100 Argon GC provides a top class service to you at all times.

Features:

Sensitivity to low ppb levels

High Frequency Argon Discharge Detector (HFADD)

Accuracy to ±1% of scale

Long Term Stability & Sensitivity

Fully Automatic

Minimum Operator involvement

Guaranteed Application



High Frequency Argon Discharge Detector (HFADD)

Application:

Analysis of Impurities in Argon

H₂, O₂, N₂, CH₄, CO & CO₂

Target Market:

- Air Separation Units
- Industrial Gas Producers
- N5, N6 & N7 Grade Argon Production
- Aerospace & Aviation Industry
- Semiconductor Industry
- Nuclear Industry

Levels of Detection

H ₂	< 15 ppb
O ₂	< 75 ppb
N ₂	< 100 ppb
CH ₄	< 30 ppb
CO	< 200 ppb
CO ₂	< 100 ppb

AGC
INSTRUMENTS

Formerly known as
GOW-MAC Instrument Co. (Ireland) Ltd

Gas Chromatography since 1965

Our Gas Chromatograph:


The AGC Series 100 HFADD GC is available in two configurations, depending on the impurities being measured. Using only components of the highest quality, enables AGC Instruments to produce a sensitive Gas Chromatograph specific for the Analysis of impurities in an Argon Matrix Gas. All fitting used are 1/8" VCR compression that ensures a perfect seal each time. The valves used are mass spec leak tested and housed in an Argon purged enclosure which provides an atmosphere free from air leaks. Thereby ensuring a contamination free environment that will provide excellent stability, sensitivity and a long working life.

Our Detector:

The gas passing through the Detector is subjected to a high electromagnetic field, thereby causing the gas to illuminate. The resulting luminance is monitored by a light dependant resistor (LDR). As the sample component elutes from the column the light intensity is altered. The output from the LDR is converted to a mV signal and is measured on the Data Capture System.

Our After Sale Support:

Installation and Commissioning support is essential to providing a complete service. Through our dedicated Distribution Network we are able to provide local support and back-up, in both native language and in English. All our Distribution partners have been selected having the technical skills required to provide an excellent support. Ongoing training is provided with specific training focused on Industrial and Specialty Gas applications.



Specifications:

Detector

Name: High Frequency Argon Discharge Detector (HFADD)
Linearity: $> 10^3$
Temp Range: 30-45°C
Sensitivity: < 100 ppb of N₂
Typical Range: 0-10ppm
Response Time: < 0.5 seconds

Power Requirements

110/220/240V, 50/60 Hz, noise and spike free

Dimensions

W = 19" ; D = 400mm ; H = 8U (360mm)

Operating Conditions

Temperature Range: +10°C to 40°C

Output Signal

0-1 V output (std) 0-10V (optional)

Gases Required

Carrier Gas: Ultra Pure Ar N6.0 ; 35-80ml/min
Actuator Gas: Clean dry air @ 3 bar

Please ask about our Argon Carrier Gas Purifier Option
Note: Special Dual Changeover Cylinder System is advised

Gas Connections

1/8" stainless steel with VCR compression fittings with 1/16" stainless steel tubing

Contact Details:

AGC Instruments Ltd, Bay K14a, Industrial Estate, Shannon, Co. Clare, Ireland

T: +353 61 471632 F: +353 61 471042 E: sales@agc-instruments.com W: www.agc-instruments.com