

## Műanyag laboratóriumi eszközök



### Spritzflakon, PE-LD



Kód	Térfogat (ml)	GL	Ø (mm)
400094588	50	18	37
400094688	100	18	43
400094988	250	25	58
400095088	500	25	74
400095188	1000	32	92

### Cseppentőflakon, PE-LD



Kód	Térfogat (ml)	GL	Ø (mm)
400094587	50	18	37
400094687	100	18	43
400094987	250	25	58
400095087	500	25	74
400095187	1000	32	92

### Erlenmeyer lombik, PP, bõnyakú, csavaros kupakkal, dugóval is használható



Kód	Térfogat (ml)	Osztás (ml)	GL	NS
400666941	50	10	40	34/35
400667941	100	20	40	34/35
400668941	250	50	52	45/40
400669941	500	100	52	45/40
400670941	1000	200	52	45/40

**Erlenmeyer lombik, PMP, bőnyakú, csavaros kupakkal, dugóval is használható**



Kód	Térfogat (ml)	Osztás (ml)	GL	NS
400066695	50	10	40	34/35
400066795	100	20	40	34/35
400066895	250	50	52	45/40
400066995	500	100	52	45/40
400067095	1000	200	52	45/40

**Folyadékflakon, PE-LD, szűknyakú, kerek forma, csavaros kupakkal**



Kód	Térfogat (ml)	GL	magasság (mm)	Ø (mm)
400094389	10	11	42	22
400094489	20	11	59	24
400094589	50	18	73	37
400094689	100	18	95	43
400094989	250	25	135	58
400095089	500	25	180	74
400095189	1000	32	215	92

**Porflakon, PE-LD, bőnyakú, kerek forma, csavaros kupakkal**



Kód	Térfogat (ml)	GL	magasság (mm)	Ø (mm)
400093389	50	32	87	39
400093489	100	32	94	47
400093789	250	45	133	58
400093989	500	45	151	76
400094189	1000	63	196	91

**Porflakon, PE-LD, bőnyakú, szögletes forma, csavaros kupakkal**



Kód	Térfogat (ml)	GL	magasság (mm)	méret (mm)
<b>400092489</b>	100	32	82	41x41
<b>400092689</b>	250	45	110	55x55
<b>400092789</b>	500	45	133	70x70
<b>400092889</b>	1000	63	166	87x87

**Folyadékflakon, PP, szűknyakú, kerek forma, kúpos felsőrész csavaros**



Kód	Térfogat (ml)	GL	magasság (mm)	Ø (mm)
<b>400100389</b>	100	18	100	52
<b>400100489</b>	250	25	132	70
<b>400100589</b>	500	25	165	87
<b>400100689</b>	1000	32	202	108
<b>400100789</b>	2000	32	245	131
<b>400100889</b>	5000	45	315	178
<b>400100989</b>	10000	63	394	222

**Porflakon, PP, bőnyakú, kerek forma, kúpos felsőrész csavaros kupakkal**



Kód	Térfogat (ml)	GL	magasság (mm)	Ø (mm)
<b>400101589</b>	100	32	96	55
<b>400101689</b>	250	45	132	73
<b>400101789</b>	500	45	172	87
<b>400101889</b>	1000	63	204	108
<b>400101989</b>	2000	63	243	131

**Folyadékflakon, PP, szűknyakú, kerek forma, kúpos felsőrész NS dugóval**



Kód	Térfogat (ml)	NS	Ø (mm)	magasság (mm)
<b>400100394</b>	100	14/23	52	106
<b>400100494</b>	250	19/26	70	138
<b>400100594</b>	500	24/29	87	172
<b>400100694</b>	1000	29/32	108	213
<b>400100794</b>	2000	29/32	131	255
<b>400100894</b>	5000	45/40	178	325

**Porflakon, PP, bőnyakú, kerek forma, kúpos felsőrész NS dugóval**



Kód	Térfogat (ml)	NS	Ø (mm)	magasság (mm)
<b>400101594</b>	100	29/32	55	111
<b>400101694</b>	250	34/35	173	144
<b>400101794</b>	500	45/40	87	183
<b>400101894</b>	1000	60/46	108	214
<b>400101994</b>	2000	60/46	131	263

**Folyadékballon, PE-LD, csappal, csavaros kupakkal, fogantyúval**



Kód	Térfogat (l)	nyak Ø (mm)	Ø (mm)	magasság (mm)
<b>400081660</b>	5	45	165	335
<b>400081662</b>	10	55	210	415
<b>400081664</b>	25	55	280	525
<b>400081666</b>	50	55	365	605

### Mérlegedény, PP, „csiszolt tetővel”



Kód	Térfogat (ml)	magasság (mm)	Ø (mm)
400080340	30	50	30
400080342	25	30	40
400080343	65	60	40
400080345	50	30	50
400080346	65	35	60
400080347	200	90	60
400080348	400	120	70

### Főzőpohár, PP, kék beosztással



Kód	Térfogat (ml)	osztás (ml)	magasság (mm)	Ø (mm)
400605081	10	2	35	25
400606081	25	5	47	31
400607081	50	10	60	40
400608081	100	20	70	49
400609081	150	20	80	56
400610081	250	50	94	68
400611081	400	50	109	77
400618081	500	100	122	88
400612081	600	100	125	91
400613081	800	100	136	98
400614081	1000	100	149	102
400615081	2000	200	183	133
400616081	3000	200	214	174
400617081	5000	500	248	185

### Főzőpohár, PMP, piros beosztással



Kód	Térfogat (ml)	osztás (ml)	magasság (mm)	Ø (mm)
400060503	10	2	35	25
400060603	25	5	47	31
400060703	50	10	60	40
400060803	100	20	70	49
400060903	150	20	80	56
400061003	250	50	94	68
400061103	400	50	109	77
400061803	500	100	122	88
400061203	600	100	125	91
400061303	800	100	136	98
400061403	1000	100	149	102
400061503	2000	200	183	133
400061603	3000	200	214	174
400061703	5000	500	248	185

**Mérőpohár, PP, domborított beosztással**



Kód	Térfogat (ml)	osztás (ml)	magasság (mm)	Ø (mm)
400446941	50	2	70	40
400447941	100	2	80	50
400440941	250	5	120	70
400441941	500	10	133	91
400442941	1000	10	170	116
400443941	2000	20	215	150
400444941	3000	50	242	170
400445941	5000	100	270	210

**Mérőpohár, PP, domborított kék beosztással**



Kód	Térfogat (ml)	osztás (ml)	magasság (mm)	Ø (mm)
400446081	50	2	70	40
400447081	100	2	80	50
400440081	250	5	120	70
400441081	500	10	133	91
400442081	1000	10	170	116
400443081	2000	20	215	150
400444081	3000	50	242	170
400445081	5000	100	270	210

## Büchner-tölcsér PP



Kód	Filter Ø (mm)	hossz (mm)	Lyuk Ø (mm)
400080437	45	95	1,0
400080438	55	113	1,0
400080439	70	145	1,5
400080440	80	165	1,5
400080441	90	180	1,5
400080442	110	210	2,0
400080443	160	280	3,0
400080445	240	350	3,0

## Tölcsértartó, PP

50-120 mm átmérőjű tölcsekhez, 8-14 mm rúdátmérőjű állványokhoz



Kód	férőhely
400080268	1
400080269	2

## Szeparálótölcsér-tartó, PP

125-500 ml-es szeparálótölcsérekhez, 8-14 mm rúdátmérőjű állványokhoz



Kód	férőhely
400080970	1

### Oliva, PE-HD, 2 részes



Kód	Olivaméret (mm)
<b>400080434</b>	3-5
<b>400080435</b>	5-7
<b>400080436</b>	7-9
<b>400080535</b>	9-12
<b>400080536</b>	11-14
<b>400080537</b>	13-6

### Oliva, PP, 2 féle méretű csatlakozóval




Kód	Olivaméret (mm)
<b>400080877</b>	4-8 / 8-12
<b>400080878</b>	4-8 / 12-16
<b>400080879</b>	8-12 / 12-16

### Olivás könyökcsatlakozó, PP



Kód	Olivaméret (mm)
<b>400081250</b>	3-4
<b>400081251</b>	4-5
<b>400081252</b>	6-7
<b>400081253</b>	8-9
<b>400081254</b>	10-11
<b>400081255</b>	12-13

### Olivás elágazó csőcsatlakozó, T forma, PP



Kód	Olivaméret (mm)
<b>400080459</b>	3
<b>400080460</b>	4-5
<b>400080461</b>	6-7
<b>400080462</b>	8-9
<b>400080463</b>	10-11
<b>400080520</b>	12-13
<b>400080521</b>	14-15

### Olivás elágazó csőcsatlakozó, Y forma, PP



Kód	Olivaméret (mm)
<b>400080464</b>	3
<b>400080465</b>	4-5
<b>400080466</b>	6-7
<b>400080467</b>	8-9
<b>400080468</b>	10-11
<b>400080525</b>	12-13
<b>400080526</b>	14-15

### Gömblobbik-tartó PP



Kód	Ø (mm)
<b>400080271</b>	160

## Pipetta állvány PP



Kód	Ø (mm)
<b>400079194</b>	230