

80-2b Electric Centrifuge, Tabletop Medical Centrifuge

Product Description

Product Description:

This machine is made from high strength plastic. It can be operated conveniently. It can be used widely to make qualitative analysis to serum, plasma and so on in hospitals and biochemistry labs. CE, ISO13485/9001/14001, SGS audited factory. Low speed table centrifuge machine features in stable operation, low noise, delicate appearance, safety-regulative balance and high precision control and it is widely used in biochemical, medical research and material units.

Feature:

1. Easy to operate, with large capacity, steady performance. Applicable for blood serum and plasma separation in laboratories.
2. As a professional manufacturer, we can offer most reasonable prices
3. More than over 20 year's history that we will provide products with high quality
4. Reliable and experienced delivery agent
5. We can develop according to customer's samples or drawings..



Specifications:

Power Supply	AC 220V/110V 50Hz/60Hz
Capacity	20ml*12/10ml*8/10ml*12
Max rcf	2325Xg
Timer	0-60min
Max speed	4000rpm
Speed indication	Pointer indication
Safety control	Auto breaking system during running

Low noise, large capacity, steady performance, widely applicable for labs, biology, chemistry, physical, hospital, food and environmental protection fields.